



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **WEREWOLF COBLATION SYSTEM SMITH & NEPHEW coablation device (Ex-Display)**

New (ex-display),  
Actual product photos,  
Made in the U.S.A. (2019),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
WEREWOLF COBLATION system  
allows functions to be controlled using the handpiece, foot controller, or touchscreen interface,  
FLOW IQ technology  
- automatically balances saline flow using coablation to optimize the flow performance of the handpiece for all tissue types,  
Controlled plasma field allows for precise removal of soft tissue in the range (100-200 µm),  
Safety  
- technology enables accurate real-time monitoring of joint fluid temperature,  
Features:  
Coagulation is available for hemostasis,  
VAC mode is a unique feature that rinses the joint without changing instruments and removes loose debris,  
HI MODE is recommended for rapid, Mass debridement of less dense, more swollen tissue.  
MEDIUM MODE is recommended for tissues of medium density, such as meniscus or healthy labrum.  
LOW MODE is recommended for tissues of higher density, such as joint cartilage.  
Consumables can be retrofitted for an additional fee.  
Includes:  
Power cord,  
Wired foot controller,  
Dimensions: 44 x 50 x 15 cm,  
Weight: 12 kg,  
Valid inspection and ready to use.  
Technical Passport (Service Report) valid for 12 months.  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market.  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market.  
Financing options (Poland only): Installments, Leasing, Loan.

