



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Medical assistant TECHMED D-03-66 (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Poland (2017),

TECHMED mobile medical assistant with removable trays and a stainless steel top shelf,

Technical specifications:

4 easily removable trays,

Stainless steel top shelf,

Side basket for accessories,

Side handle for easy transport,

Mobile assistant placed on 4 wheels (2 with brakes),

Included:

4 trays

Dimensions:

Top shelf: 60 x 36 cm,

Tray: 28 x 35 cm,

Total dimensions: 98 x 44 x 88 cm,

Weight: 15 kg,

Valid inspection and ready to use,

Issued Technical Passport (Service Report) valid for 12-month warranty

Warranty:

6 months for the domestic market (Poland),

3 months for the international market

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets

Financing options (Poland only): Installments, Leasing, Loan

