



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



MEDEN-INMED SOLMED UNO therapeutic irradiation lamp (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Poland,
Power supply: 230 V,
Frequency: 50/60 Hz,
Bulb: 375 W,
Solmed UNO

is a lamp designed for performing skin tissue heating treatments using infrared thermal energy emitted by a single bulb. Additionally, depending on the desired treatment effect, a red or blue color filter can be used, which changes the radiation characteristics. It is equipped with a single tube. The lamp is used to treat certain chronic and subacute inflammatory conditions and pain syndromes.

Control panel:

Display -

located in the central part of the control panel. It displays the lamp's operating parameters.

Buttons for changing operating parameters (bulb power and operating time) - on the left side of the panel are buttons for setting the bulb power, and on the right side are buttons for setting the treatment time.

START and STOP buttons -

The START button is used to start and stop the treatment (PAUSE mode), and the STOP button is used to stop the treatment. The START button is also used to exit sleep mode.

PROGRAM button -

This button is used to enter treatment sequence programming mode. Programming mode is indicated by a LED located under the button.

Specifications:

Treatment timer: 1 - 30 min.,

Brightness adjustment: 10-100%, 10% increments,

Applications:

Inflammation of cranial nerves,

Inflammation of joints,

Post-inflammation of intercostal nerves,

Post-burn conditions,

Periarthritis,

Maxillary sinusitis,

Contractures, scars, adhesions,

Rosacea,

Adolescent acne,

Shingles,

Frostbite,

Paresthesia,

Myalgia,

Difficult-to-heal wounds,

Includes:

Red filter,

Blue filter,

Power cord,

Green and orange filters can be added for an additional fee.

Dimensions:

Overall: 50 x 67 x 172 cm,

Weight: 14.5 kg,

It has a current inspection and is ready to use,

A Technical Passport (Service Report) is issued, valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,

The warranty period does not cover light bulbs and other light sources installed in the device,

Financing options (Poland only): Installments, Leasing, Loan,

