



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



JORDAN H gynecological and urological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the Czech Republic,
Medically certified gray upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
JORDAN H gynecological and urological chair - designed for gynecological and urological examinations and procedures in healthcare facilities,
Construction made of powder-coated steel profiles,
Specifications:
Manual adjustment:
Backrest: 48°,
Seat: 15°,
Bilateral armrests with freely adjustable spacing and height,
Sliding waste bowl with dimensions: 40 x 18 x 15 cm,
Paper holder with a width of 52 cm,
Retractable patient step,
Includes:
2 pcs - leg rests Patient,
Dimensions:
Backrest: 60 x 78 cm,
Seat: 60 x 40 cm,
Total dimensions: 80 x 115 x 142 cm,
Weight: 50 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

