



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **JORDAN H gynecological and urological chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in the Czech Republic,  
Medically certified gray upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),  
JORDAN H gynecological and urological chair - designed for gynecological and urological examinations and procedures in healthcare facilities,  
Construction made of powder-coated steel profiles,  
Specifications:  
Manual adjustment:  
Backrest: 48°,  
Seat: 15°,  
Bilateral armrests with freely adjustable spacing and height,  
Sliding waste bowl with dimensions: 40 x 18 x 15 cm,  
Paper holder with a width of 52 cm,  
Retractable patient step,  
Includes:  
2 pcs - leg rests Patient,  
Dimensions:  
Backrest: 60 x 78 cm,  
Seat: 60 x 40 cm,  
Total dimensions: 80 x 115 x 142 cm,  
Weight: 50 kg,  
Updated inspection and ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

