



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



UMF MEDICAL EXAM aesthetic medicine chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Upholstery with medical certificate in blue color (it is possible to replace the upholstery with a new one, from our color pattern for an additional fee),
UMF MEDICAL EXAM chair
is a multifunctional chair ideally suited for ENT and ophthalmology surgeries and can also serve as a blood collection chair,
Through its built-in features, it is also ideal for aesthetic medicine and for performing various types of procedures,
Adjustment of the table's functions is carried out with a handheld remote control,
Technical specifications:
Electric adjustment by handheld remote control:
Backrest adjustment in the range: 52° (simultaneous synchronization of the leg rest with the patient's armrests),
Height adjustment in the range: 61 - 91 cm,
Manual retractable patient armrests on both sides,
Retractable stainless steel tray for accessories with dimensions: 20.5 x 30 cm,
Includes:
Wired handheld remote control,
Power cord,
Dimensions:
Back rest: 42 x 85 cm,
Seat: 46 x 47 cm,
Leg rests: 43 x 29 cm,
Armrests: 11.5 x 58 cm,
Base: 52 x 97 cm,
Total: 80 x 97 x 140/170 cm,
Maximum load capacity: 170 KG,
Weight: 100 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

