



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER VAPOR 2000 (Sevoflurane/Sevoflurane) inhalant vaporizer (Reconditioned)

Product temporarily unavailable

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany (2017),

Drager Vapor 2000

is a safe and precise vaporizer for inhalation agents.

Product features:

The vaporizer's volume is 300 ml, allowing the entire contents of a standard 250 ml anesthetic bottle to be poured into the vaporizer, even if anesthetic remains inside from the previous refill.

The vaporizer does not require calibration during use.

Unique transport feature - vaporizers can be carried in any position, as the transport lock prevents anesthetic spillage, so they can be used immediately after transport.

All filling systems have been designed specifically for specific inhalation anesthetics, eliminating the risk of mistakenly using the wrong anesthetic agent.

In a given The vaporizer features the Drager DW-2000 mount.

Infusion type: Drager Fill.

Sevoflurane concentration can be adjusted from 0 to 8%.

Includes:

Instructions for use in English (PDF).

Dimensions: 11 x 25 x 20 cm.

Weight: 9 kg.

Cat 03.

Updated inspection and ready for use.

Issued Technical Passport (Service Report) valid for 12 months.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

