



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



MEDIFA MOT 6000 operating table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production, (2013),
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 15 Ah,
- Frequency: 50/60 Hz,
- Power: 300 VA,
- 6-segment,
- **Upholstery with medical certificate in green,**
- **MEDIFA MOT 6000** is a comfortable, mobile operating table with electrohydraulic drive, designed for all surgical specialties,
- **Electrically adjustable using a wired remote control and a control panel on the table column in many planes:**
 - Height range: 69 – 105 cm,
 - Trendelenburg and anti-Trendelenburg positions in the range of: 30°,
 - Lateral tilt in the range: 20°,
 - Lower backrest angle range: from -35° to 70°,
 - Backrest angle of the upper back section in the range: from -25° to 55°,
 - Longitudinal travel range: 31 cm,
 - Wheel lock,
 - Position "0",
- **Manual adjustment:**
 - Headrest inclination angle range: from -35° to 55°,
 - Split leg segment, sideways, adjustable inclination angle,
- It has X-RAY transparency,
- Suitable for C-arm,
- Rails for mounting additional equipment,
- Colored levers for quick and easy table modification,
- The base is made of stainless steel, profiled for good access for the operators' feet by narrowing the housing downwards,
- Mobile on four double wheels with a diameter of 12 cm, with the option of directional driving and central locking using a remote control,
- **Included in the set:**
 - Power cord,
 - Wired remote control,
- Dimensions:
 - Bed base: 210 x 55,
 - Overall: 210 x 60 x 69/105 cm
- Maximum load capacity: 300 KG,
- Weight: 190 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

