



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS PAGEWRITER TRIM III ECG device (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the U.S.A.
Power supply: 230 V,
Battery power: via a new built-in 12 V / 2.3 Ah battery,
Frequency: 50/60 Hz,
PHILIPS PAGEWRITER TRIM III
is a compact and ergonomic cardiograph thanks to its simplified functionality.
The cardiograph has a high-resolution color display.
Features:
12-lead ECG (adults and pediatrics),
PHILIPS' 12-lead algorithm is used for interpretive analysis of waveform amplitudes, durations, and morphology, as well as the associated rhythm, for adult and pediatric patients.
Automatic and configurable pacemaker detection.
Real-time ECG signal quality indicator.
On-screen ECG report preview.
Specifications:
ECG 12-lead,
Alphanumeric keyboard,
Display: Color TFT 640 x 480 px,
Sampling rate: 2000 samples per second per electrode/lead,
Built-in memory for 50 studies,
Built-in PCMCIA card reader,
Report formats:
3x4 1R, 3R (Standard, Cabrera),
6x2 (Standard, Cabrera),
Panoramic 12 (Cabrera),
12x1 (Standard, Cabrera),
Rhythm (up to 12 selected leads),
12-wire extended measurements,
(1 minute of continuous waveform data for 1 selected lead),
Automatic frequency response:
0.05-150 Hz / 0.15-150 Hz / 0.5-100 Hz / 0.05-100 Hz / 0.15-100 Hz / 0.5-100 Hz / 0.05-40 Hz / 0.15-40 / 0.5-40 Hz,
Rhythmic frequency response:
0.05-150 Hz / 0.15-150 Hz / 0.05-150 Hz / 0.15-100 Hz / 0.05-40 Hz / 0.15-40 Hz,
Built-in thermal printer:
4-channel, real-time printout (paper width: 210 mm x 300 mm),
Menu language: English,
Includes:
Power cord,
10-lead ECG cable,
6 EPP suction cup electrodes,
4 EKK clip limb electrodes,
Weight: 7.3 kg,
Cat. 02,
It has a current inspection and is ready for operation.
A Technical Passport (Service Report) is issued and valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
The warranty can be extended for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan.

