



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STIEGELMEYER P-BETT ELVIDO rehabilitation bed (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany (2020),

Power supply: 230 V,

Frequency: 50/60 Hz,

The Stiegmeyer Elvido P-Bett

is a modern electric rehabilitation bed designed with patient comfort and caregiver ergonomics in mind. It also boasts an elegant appearance, durable construction, and extensive functionality.

The bed is designed to provide patients with maximum safety and comfort in everyday use. Electric adjustment allows the bed to be adjusted to multiple positions – from resting to rehabilitation – which supports the treatment and care process.

Thanks to its very low starting height (approx. 25 cm), the bed minimizes the risk of fall-related injuries, which is particularly important when caring for seniors or patients with dementia. The modern Vario Safe side rail system allows for quick adjustment of their height and length to the user's needs.

Control is via a handheld remote control with the option to lock individual functions, increasing safety.

Technical specifications:

Height adjustment range: 25 – 80 cm,

Backrest angle range: 70°,

Thigh section angle: up to 35°

Functions: Trendelenburg and Reverse Trendelenburg positions, seat position, electric adjustment of the four-segment bed surface,

Mobility: Wheels with a central locking system for easy transport,

Construction: painted steel Powder-coated, decorative wooden panels.

Includes:

Patient arm with handle,

Wired remote control,

Power cord,

Maximum load capacity:

Working: 225 kg,

Maximum patient weight: 185 kg,

Dimensions:

Sleeping area: 90 x 200 cm,

Weight: 150 kg,

Updated inspection and ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan.

