



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHYSIO-CONTROL docking station (charger) for LIFEPAK 15 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Canadian production,
- Power supply: 230 V and via car connector 13.2 V,
- Frequency: 50/60 Hz,
- **Designed for charging Lifepak 15 lithium-ion batteries,**
- **Two battery slots,**
- **Included in the set:**
 - Power cord,
 - 13.2V car charger,
- Dimensions: 18.5 x 27.5 x 11 cm,
- Weight: 3 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

