



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diagnostic and surgical table OOSTWOUD (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Dutch production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Oostwoud 2-segment diagnostic and examination table,
Medical approved upholstery in black color after purchase will be replaced with a new one in the color of customer's choice (from our color chart),
Specification:
Frame made of powder-coated profiles,
Electric height adjustment with wired remote control in the range: 54-90 cm,
Manual adjustment of the angle of the headrest in the range: 0 - 55°,
Paper pad hanger with a width of 55 cm,
Mobile on 15 cm diameter casters with central brake and directional travel option,
Included:
Power cord,
Wired remote control,
Dimensions:
Recliner: 195 x 75 cm,
Overall: 195 x 75 x 54/90 cm,
Weight: 85 kg,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

