



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Synoptophore OCULUS 58100 (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany,

Power supply: 230 V,

Frequency: 50/60 Hz,

Technical specifications:

Light intensity adjustment: From 10% to 100% (adjustable in 5% steps),

Image adjustment range on the vertical axis: +/- 20°,

Examination distance adjustment range (PD): 42 - 73.5 mm,

Flash adjustment: automatic / manual,

Menu language: German,

Included:

Power cord,

20 pairs of slides,

Overall dimensions: 60 x 30 x 40 cm,

Weight: 25 kg,

Updated inspection and ready to use,

Issued Technical Passport (Service Report) valid for 12-month warranty

Warranty:

6 months for the domestic market (Poland),

3 months for the international market

Possibility of extending the warranty to 12 months for the domestic and international markets for an additional fee

Warranty period does not cover light bulbs and other light sources installed in the device

Financing options (Poland only): Installments, Leasing, Loan

