



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NIDEK NT-2000 non-contact tonometer (MIII) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Japan,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 70 VA,
Technical specifications:
Measuring range: 1-60 mmHg,
Built-in display: 5" TFT LCD,
Working distance: 11 mm (from the air nozzle to the cornea),
Electric bowl lifting,
Automatic safety lock to prevent contact with the patient's eye,
Computer interface for data transmission,
Manual chin rest adjustment within a range of 7 cm,
Menu language: English,
Includes:
Instruction manual in Polish and English,
Power cord,
Dimensions: 26 x 48.7 x 45.7 cm,
Weight: 17 kg,
Includes a current Inspected and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
The warranty period does not cover light bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

