



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HETTICH MIKRO 200 laboratory centrifuge (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 240 VA,
Features:
Metal housing and lid,
Inspection window in the lid,
One-handed lid lock,
Lid fall protection,
Emergency lid lock release,
Stainless steel chamber,
Automatic rotor recognition,
Brushless drive,
Error display,
Temperature display in °C and °F,
Switch-off switch in case of imbalance,
Illuminated panel with current values of all parameters,
Audible message at the end of the spin cycle,
Specifications:
RPM (revolutions per minute): From 500 to 15,000 (adjustable in increments of 10),
Max. RCF: 21.382,
Pulse button for short cycle mode,
Easy operation via keypad and knob,
4 program memories for greater customization,
9 individual acceleration and deceleration stages,
Menu language: German,
Included:
Operating instructions in ENG, GER, FR, IT (PDF),
Power cord,
HETTICH 2424 rotor (24 x 3g),
Brand new HETTICH rotor cover,
Dimensions: 27.5 x 34.4 x 26 cm,
Weight: 11.5 kg,
Updated inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee. 12 months for domestic and international markets. Financing options (Poland only):
Installments, Leasing, Loan.

