



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



CARDIETTE AR600 ECG device (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Italian production,
Power supply: 230 V,
Battery backup: via a new built-in 1.2 V / 2000 mAh battery,
Frequency: 50/60 Hz,
The AR600 ECG is a portable electrograph with a dual power supply.
Features:
ECG recording in automatic or manual mode,
ECG printing on single-channel 60 mm wide paper using high-resolution thermal printing,
Save the latest ECG recording in automatic mode and print additional copies of the ECG as needed,
Technical specifications:
ECG signal: 12 leads,
Recording sensitivity:
Manual: 5 / 10 / 20 mm/mV +/- 5% ,
Automatic (depending on the number of printed channels),
Recording speed: 25 – 50 mm/s,
Includes:
User manual in English (PDF),
Power cord with power supply,
10-lead ECG cable,
Dedicated trolley,
Dimensions: 25 x 18.5 x 6 cm,
Weight: 1 kg,
Updated inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

