



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEISEGANG 3MLF Colposcope (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Halogen bulb: 15 V / 150 W,
Technical specifications:
Beam path: Convergent,
Eyepieces: direct viewing,
Magnification: 7.5x / 15x / 30x,
Field of view diameter: 46 / 23 / 11.5 mm,
Light field diameter: 58 mm,
Diopter correction: +7 to -7 diopters,
Horizontal fine adjustment: 40 mm,
Vertical adjustment: 80 mm,
The colposcope is equipped with a light source with the ability to freely adjust the intensity using knobs,
Colposcope mounted on a 2-segment frame with free rotation adjustment,
Includes:
Power cord,
Dimensions:
Frame height: 94 cm,
Weight: 18 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
The warranty period does not cover light bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,



