



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRECISE PS-800A Smoke Evacuator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Manufactured by P.R.C. (2018),
Power supply: 230 V,
Frequency: 50/60 Hz,
The PRECISE PS-800A

is a smoke evacuation system that captures particles and adsorbs gases from surgical smoke. It has been specifically designed to improve visibility and reduce potential health risks associated with surgical smoke.

It can be used in both open and laparoscopic procedures.

Flow rate can be adjusted in 5 steps.

Features:

Designed with a centrifugal pump with high suction power and flow.

Easy to set up and use.

LED display showing filter life and suction power.

Simple filter replacement (the filter can capture particulates and microorganisms down to 0.012 µm with 99.99% efficiency).

Used with a CO2 laser system during surgery.

Equipped with LED indicators for:

Filter status,

Flow rate,

Service required.

Includes:

New filter,

Foot switch,

Suction tube,

Power cord.

Dimensions: 18 x 44 x 30 cm.

Weight: 10 kg,

It has a current inspection and is ready for operation.

A Technical Passport (Service Report) is issued, valid for 12 months.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

The warranty can be extended for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

