



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BKG BA-4 contrast meter (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany,

Power supply: 230 V,

Frequency: 50/60 Hz,

The BA-4 contrast meter is a pre-programmed device designed to test contrast vision and glare sensitivity.

Tests can be performed separately for each eye; eyes can be tested simultaneously or individually.

Includes:

Instruction manual in Polish,

Power cord with power supply,

Dimensions: 27 x 40 x 40 cm,

Weight: 2 kg,

It has a current inspection and is ready to use.

Issued Technical Passport (Service Report) valid for 12 months.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

Possibility to extend the warranty for an additional fee. 12 months for domestic and international markets.

The warranty period does not cover light bulbs and other light sources installed in the device.

Financing options (Poland only): Installments, Leasing, Loan.

