



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PROMA REHA PLE-N85-T rehabilitation bed (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in the Czech Republic,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 350 W,  
PROMA REHA PLE-N85-T  
is a comfortable and mobile rehabilitation bed with electric adjustment.  
Specifications:  
Electric adjustment via a wired remote control:  
Height adjustment range: 42 – 82 cm,  
Backrest adjustment range: 78°,  
Legrest adjustment range: 52°,  
Retractable side rails on both sides,  
Rubber bumpers on each side of the bed,  
Built-in emergency CPR safety system,  
Central foot brake,  
Mattress not included – mattress can be added for an additional fee,  
Included:  
Wired remote control Manual,  
Dimensions:  
Bedside table dimensions: 72 x 190 cm,  
Backrest dimensions: 72 x 88 cm,  
2-segment legrest: 72 x 70 cm,  
Total dimensions: 82 x 210 cm,  
Weight: 90 kg,  
Updated inspection and ready for use,  
Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

