



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynecological chair MAQUET RADIUS 1557.02AJ (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 290 VA,
Medically certified blue upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
Electromechanical adjustment via wired hand and foot controls:
Bed height within a range of 61 - 91 cm,
Tilt angle:
Back section within a range of (- 9°) / 0° / 57°,
Seat within a range of 18°,
When adjusting the back section angle, the seat angle automatically adjusts,
3 pre-programmed buttons on the hand remote control,
S button to stop the programmed mode button,
Electric reclining of the seat into a lying position,
The seat features:
A 38 cm wide roll storage compartment for a roll of backing paper,
A cutout in the back section upholstery for backing paper,
Retractable footrests,
The seat can be reclined into a lying position,
Includes:
Instruction manual in Polish (PDF),
Power cord,
Wired handheld remote control,
Wired foot pedal remote control,
Waste bowl,
Option to retrofit an additional colposcope mount for an additional fee,
Dimensions:
Back section: 91 x 61 cm,
Seat area: 61 x 31 cm,
Length in the unfolded position: 175 cm,
Overall dimensions: 155 x 76 x 138/168 cm,
Maximum load capacity: 180 kg,
Weight: 94 KG,
Cat 02,
It has a current inspection and is ready for operation.
A Technical Passport (Service Report) is issued, valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan.

