



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RENALCHAIR CE 300 dialysis chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 160 VA,
Medically certified upholstery in blue and navy blue (upholstery can be replaced with a new one from our color chart for an additional fee),
The RENALCHAIR CE 300 dialysis chair
provides flexible adjustments with a wide range of positioning options.
The chair offers excellent mobility, making it suitable for easy movement around medical facilities.
Technical specifications:
Adjustment via wired handheld remote control:
Seat height adjustment range: 57.5 – 71 cm,
Backrest tilt adjustment range: 87° / 0°,
Free adjustment Legrest tilt range: -40° / 0° / 25°,
Trendelenburg position possible,
Synchronized backrest and legrest adjustment possible,
2 built-in pre-programmed positioning buttons,
Bilateral armrests with manual free rotation,
Manually adjustable patient step extension,
Built-in paper holder,
Mobile chair with built-in central brake,
The chair upholstery shows visible signs of wear that do not affect comfort,
Includes:
Power cord,
Patient headrest,
Dimensions:
Backrest: 54 x 85 cm,
Seat: 54 x 46 cm,
Legrest: 54 x 61 cm,
Armrests: 16 x 65 cm,
Maximum load capacity: 180 kg,
Weight: 80 kg,
Current inspection and It is ready for operation.
A Technical Passport (Service Report) is issued and is valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
The warranty can be extended for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan.

