



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **RENALCHAIR CE 300 ENT chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 160 VA,  
Medically certified upholstery in blue and navy blue (upholstery can be replaced with a new one from our color chart for an additional fee),  
The RENALCHAIR CE 300 ENT chair  
ensures easy adjustment with a wide range of positioning options.  
The chair offers excellent mobility, making it suitable for easy movement around medical facilities.  
Technical specifications:  
Adjustment via wired handheld remote control:  
Seat height adjustment range: 57.5 – 71 cm,  
Back tilt adjustment range: 87° / 0°,  
Free adjustment Legrest tilt range: -40° / 0° / 25°,  
Trendelenburg position possible,  
Synchronized backrest and legrest adjustment possible,  
2 built-in pre-programmed positioning buttons,  
Bilateral armrests with manual free rotation,  
Manually adjustable patient step extension,  
Built-in paper holder,  
Mobile chair with built-in central brake,  
The chair upholstery shows visible signs of wear that do not affect comfort,  
Includes:  
Power cord,  
Patient headrest,  
Dimensions:  
Backrest: 54 x 85 cm,  
Seat: 54 x 46 cm,  
Legrest: 54 x 61 cm,  
Armrests: 16 x 65 cm,  
Maximum load capacity: 180 kg,  
Weight: 80 kg,  
Current inspection and It is ready for operation.  
A Technical Passport (Service Report) is issued and is valid for 12 months.  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market.  
The warranty can be extended for an additional fee to 12 months for the domestic and international markets.  
Financing options (Poland only): Installments, Leasing, Loan.

