



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HILL-ROM AFFINITY 4 Birthing Bed (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 7.0 Ah,
- Frequency: 50/60 Hz,
- **The Affinity 4** is an electrically operated maternity bed for gynecology, obstetrics, and perinatology departments. The bed combines safety, comfort, and ease of use with unrivaled durability and reliability.
- **Technical specifications:**
 - Backrest segment adjustment: 0° / 68°,
 - Maximum seat segment lift: 0° / 15°,
 - Trendelenburg position: 8°,
 - Electric height adjustment in the range: 56 – 97 cm,
 - Central control and braking,
 - Immediate resuscitation position (CPR function),
 - Lock button and indicator,
 - Birthing handles,
 - OneStep drop-down railings,
 - Trendelenburg positioning mechanism,
 - Foot rests,
- **Medically certified mattress in blue,**
- **Possibility of retrofitting a dedicated bowl for an additional fee,**
- **Included in the set:**
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 229 x 99 cm,
- Maximum load capacity: 227 KG,
- Weight: 218 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

