



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MERIVAARA ADATTO transport trolley (Reconditioned)

Reconditioned (used),
Technical condition: good,
Visual condition: good,
Actual product photos,
Made in Finland,
Power supply: 230 V,
Frequency: 50/60 Hz,
Specifications:
Electric adjustment via wired remote control:
Height adjustment: 46 – 81 cm,
Backrest: 0° / 70°,
2-segment legrest adjustment,
Manual adjustment of Trendelenburg and Reverse Trendelenburg positions,
Retractable side rails on both sides,
Bumpers located on each side of the stroller,
Built-in central brake,
Dimensions:
Backrest: 72 x 78 cm,
Seat: 72 x 20 cm,
2-segment legrest: 72 x 84 cm,
Maximum weight capacity: 225 kg,
Weight: 80 kg,
It has a current inspection and is ready for operation.
A Technical Passport (Service Report) is issued, valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan.

