



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MRS 2000 magnetotherapy device (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
MRS 2000 Magnetic Mattress  
operates on the principle of therapeutic magnetic fields. Thanks to its unique design, it generates and introduces a pulsating resonant magnetic field at the appropriate frequency into the body, which reaches all body cells, neutralizing the harmful effects of electromagnetic smog. Activates the body's self-regulation mechanism,  
Magnetic field intensity can be adjusted to 10 levels,  
Applications:  
Stress (depression, insomnia, neurosis),  
Immune system disorders,  
Chronic respiratory infections,  
Circulatory system impairment,  
Soft tissue damage,  
Accelerates wound healing and pressure sores (diabetic wounds),  
Accelerates the healing of bone fractures in degenerative diseases,  
After injuries (limb sprains, dislocations, contusions),  
Useful in rehabilitation after strokes, hemorrhages, multiple sclerosis, atherosclerosis,  
Benefits:  
Improves peripheral blood circulation,  
Accelerates the absorption of hematomas and lymphatic edema,  
Relieves pain of rheumatic origin,  
Relieves neuralgia,  
Reduces nervous tension, fatigue, and anxiety,  
Provides Restful sleep,  
Minimizes the effects of stress,  
Reduces post-exercise restitution (increases oxygen uptake),  
Regulates blood pressure and has a positive effect on cardiovascular diseases,  
Reduces migraine headaches,  
Beneficial in the treatment of diabetes, atherosclerosis, and osteoporosis,  
Accelerates the healing of wounds, bone fractures, and joint damage,  
Reduces the deposition of cholesterol plaques in blood vessels,  
Detoxifies the body and increases its resistance to disease,  
Combats inflammation,  
Includes:  
Power cord,  
2 x applicator,  
Dimensions of two applicators:  
1: 65 x 155 cm,  
2: 50 x 28 cm,  
Total weight: 5 kg,  
Cat02,  
Updated inspection and ready to use,  
Issued Technical Passport (Service Report) valid for 12-month warranty  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets  
Financing options (Poland only): Installments, Leasing, Loan

