



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ MEDI-MATIC 115.515 gynecological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
New upholstery with medical certificate in green (upholstery can be replaced from our color chart for an additional fee),
Schmitz Medi-Matic
is a chair used for gynecological and urological examinations and procedures,
Smooth surface for optimal cleaning and hygiene,
Smooth chair adjustment using a foot switch:
Seat height range: 74 - 95 cm,
Back section inclination range: - 16° / 0° / 43°,
Pelvic section inclination range: - 51° / 0° / 33°,
The chair features:
Paper space: 40 cm wide cm,
Manually adjustable patient handles,
Rails for mounting additional accessories,
Manually adjustable Goepel knee rests,
Sliding stainless steel tray with dimensions: 29.5 x 32.5 x 6 cm,
Manually extendable thigh support,
2 built-in 230 V sockets for powering additional accessories,
Included:
Power cord,
Foot switch,
2 knee rests,
2 patient handles,
Headrest,
Waste bowl,
Dimensions:
Headrest: 30 x 20 x 12 cm,
Backrest: 58 x 87 x 5 cm,
Seat: 56 x 32 x 5 cm,
Base: 54 x 94 cm,
Overall: 122 x 90 x 133/154 cm,
Maximum weight capacity: 135 kg,
Weight: 174 kg,
It has a current inspection and is ready for operation.
A Technical Passport (Service Report) is issued, valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan

