



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **TOPCON SL-2D Slit Lamp (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Japan,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Bulb: 6 V / 20 W,  
Technical specifications:  
Eyepiece: 12.5x,  
Magnification adjustment: 10x / 16x / 25.6x,  
Smooth slit width adjustment,  
Smooth slit height adjustment,  
Slit angle adjustment: 0 - 180°,  
Rotating eyepiece (lamp): 90° right / left,  
Freely adjustable eyepiece spacing,  
Built-in filters:  
Blue,  
Green,  
Gradual light intensity adjustment,  
Manual adjustment of the patient's chin rest within the range of: 8 cm,  
Point Fixation,  
The equipment is placed on an electric MDT table with electric height adjustment in the range of: 70 – 90 cm,  
Overall dimensions: 85 x 50 x 130 cm,  
Weight: 20 kg,  
It has a current inspection and is ready for operation,  
A Technical Passport (Service Report) is issued, valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

