



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NEW HOSPIVAC 350 CA-MI FS2 surgical suction unit (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Italian production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 230 VA,
The NEW HOSPIVAC 350 FS2 suction unit
is used to remove various secretions from internal organs. Its great advantage is that it operates continuously and does not lose efficiency. Its quiet operation makes it suitable for use in the operating room or intensive care unit.
The device is equipped with a special pump that supports suction. Thanks to the two-wheel lock, the product stands firmly in place.
Technical specifications:
Maximum vacuum: -90 kPa,
Maximum flow: 60 l/min,
Noise level: 51.7 dBA,
Duty cycle: continuous operation,
Includes:
Power cord,
2x 2-liter autoclavable jars,
Wired foot switch,
Dimensions: 46 x 42 x 85 cm,
Weight: 13 kg,
Cat 03,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing available (only Poland): Installments, Leasing, Loan,

