



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybierański  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **STATIM 2000S autoclave (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Canada,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
STATIM 2000S

is a fast and effective cassette autoclave (for 2-liter cassettes). The sterilization process is fully automatic; simply select one of the available programs to effectively sterilize dental and surgical instruments, medical equipment, and dressings, among others.

STATIM 2000S solutions offer comprehensive sterilization options for instruments with or without packaging. STATIM Dri-Tec technology ensures rapid sterilization and drying.

STATIM 2000S

uses convection heating to dry instruments, utilizing the heat remaining in the system after the sterilization phase. Heat is collected and released into the cassette to rapidly dry properly loaded instruments in the STATIM cassette.

Cycles:

EMPTY UNWRAPPED instruments (S) 134°C / 3.5 min.,

EMPTY UNWRAPPED instruments (S) 134°C / 18.0 min.,

SOLD UNWRAPPED instruments (N) 134°C / 3.5 min.,

EMPTY WRAPPED instruments (S) 134°C / 3.5 min.,

EMPTY WRAPPED instruments (S) 134°C / 18.0 min.,

RUBBER & PLASTIC (S) 121°C / 15.0 min.,

RUBBER & PLASTIC (S) 121°C / 30.0 min.,

Drying: 60 min.,

Menu in Language: Polish

Includes:

Instruction manual in English (PDF),

Power cord,

STATIM sterilization cassette,

Condensate container,

Dimensions:

Overall: 41.5 x 48.5 x 15 cm,

External cassette dimensions: 19.5 x 41 x 4 cm,

Internal cassette dimensions: 28 x 18 x 4 cm,

Weight: 21 kg,

Updated inspection and ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

