



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SCHMITZ MEDI-MATIC gynecological chair 115.505 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Medically certified upholstery in white (upholstery can be replaced from our color chart for an additional fee),  
Schmitz Medi-Matic is used for gynecological and urological examinations and procedures,  
Technical specifications:  
Electromechanical adjustment via foot switch:  
Seat height: 66 - 100 cm,  
Tilt angle:  
Back section: (- 15°) / 0° / 50°,  
Possibility of reclining the chair to a lying position,  
Goepel-type knee rests, adjustable in many dimensions and attached to the armrests (2 pcs),  
Equipment:  
Double-sided rails for mounting additional Accessories,  
Manually adjustable headrest,  
Goepel knee rests adjustable in multiple dimensions and mounted on side rails - 2 pcs,  
2-function foot control,  
Extendable patient step,  
Extendable waste pan with dimensions: 32 x 30 x 6 cm,  
Included:  
Power cord,  
Wired attachment,  
Dimensions:  
Back section: 87 x 57 x 8 cm,  
Seat: 55 x 32 x 8 cm,  
Total length in the lying position: 140 cm,  
Overall dimensions: 150 x 135 x 86 cm,  
Maximum weight capacity: 135 kg,  
Weight: 173 kg,  
Up-to-date inspection and ready for use,  
Issued Technical Passport (Service Report) valid for 12-month warranty  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets  
Financing options (Poland only): Installments, Leasing, Loan

