



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **BAKMED PH121 treatment lamp (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Poland,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
BAKMED PH121 mobile lamp  
intended for use in various types of treatments,  
The lamp is placed on an adjustable 65 cm arm with adjustable head height from 77 to 152 cm,  
The 5-arm mobile lamp base facilitates comfortable use,  
Weight: 12 kg,  
It has a current inspection and is ready to operate,  
A Technical Passport (Service Report) is issued, valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,  
Period The warranty does not cover light bulbs and other light sources installed in the device.  
Financing options (Poland only): Installments, Leasing, Loan.

