



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GREINER operating table (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany (2023),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Medically certified upholstery in dark gray (upholstery can be replaced with a new one from our color chart for an additional fee),  
The GREINER operating table  
is designed for general surgery procedures and surgeries.  
The table has a simple design and, thanks to its electric adjustment using a foot control, improves comfort of use.  
It has a mobile base with a built-in central brake.  
Specifications:  
Electric adjustment using a foot control:  
Electric backrest tilt: 90° / 0° / -7°,  
Electric height adjustment in the range: 58 - 88 cm,  
Manual adjustment:  
Headrest height adjustment in the range 15 cm,  
Free headrest tilt adjustment,  
Tilt adjustment of the 2-segment footrest within a range of 90°,  
Manual table rotation within a range of 360°,  
Retractable armrests on both sides,  
Built-in central brake,  
Included:  
Wired foot control,  
Dimensions:  
Headrest: 28 x 22 cm,  
Backrest: 50 x 58 cm,  
Seat: 54 x 52 cm,  
Maximum load capacity: 200 kg,  
Weight: 95 kg,  
Updated inspection and ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for up to 12 months of subsidy for domestic and international markets. Financing options (Poland only):  
Installments, Leasing, Loan.

