



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CONSULTRONIX FOTOFINDER MEDICAM 1000S VEXIA video-dermatoscope (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Illustrative + Actual product photos,  
Made in Germany (2021),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 350 VA,  
Software: FOTOFINDER UNIVERSE 3.4.3.0,  
The VEXIA model is a versatile video dermoscopic system that allows for both quick skin scans for signs of melanoma and long-term photodocumentation of moles, lesions, skin diseases, hair, and nails.  
The device guarantees a crystal-clear image of structures at optical magnification up to 120x (standard D-scope IV lens).  
This versatile video dermoscopic system allows for both quick skin scans for signs of melanoma and long-term photodocumentation of moles, lesions, skin diseases, hair, and nails. The device guarantees a crystal-clear image of structures under optical magnification, introducing the user to a high-magnification dermatoscope.  
Thanks to its ergonomic, high-quality design, video documentation is very comfortable.  
Advantages:  
Full HD CrystalView image,  
Precise optical system,  
Built-in polarized light,  
Immersion imaging mode,  
D-Scope IV lens provides up to 140x magnification.  
Applications:  
Videodermoscopy,  
Trichoscopy,  
Onychoscopy,  
Capillaroscopy,  
Clinical – macroscopic skin images.  
The set includes:  
EIZO EV2785 monitor,  
MEDICAM 1000s camera,  
1000s docking station,  
D-Scope IV lens,  
FOTOFINDER MEDICAL SERVER computer,  
VEXIA II cart,  
Keyboard Wired,  
Wired mouse,  
Weight: 58 kg,  
Updated inspection and ready to use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,  
The warranty period does not cover light bulbs and other light sources installed in the device,  
Financing options (Poland only): Installments, Leasing, Loan,

