



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE ICC 300 surgical diathermy (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 300 W,
Smooth power regulation,
Digital power setting displays,
Easy-to-use, laminated control panel,
Control panel keyboard in EN, DE,
It has the following functions:
Monopolar cutting - max. 300 W:
AUTO CUT - cutting with automatic voltage control,
HIGH CUT - cutting tissues with poor electrical conductivity,
Four coagulation types - max. 120 W:
SOFT COAG - precise control of the coagulation process and prevention of carbonization through coagulation mode,
FORCED COAG - when fast and effective coagulation is needed,
SPRAY COAG - suitable for coagulation with argon gas,
AUTO BIPOLAR - bipolar coagulation,
Includes:
Brand new monopolar electrode cable with BOWA ERGO PEN holder - 2 buttons, 2.4 mm shaft
(electrodes not included, can be retrofitted for an additional fee),
Brand new cable with BOWA clip for disposable neutral electrodes
(electrodes not included, can be retrofitted for an additional fee),
Brand new bipolar cable, BOWA forceps
(bipolar tweezers can be retrofitted for an additional fee),
Brand new foot switch,
Instruction manual in English (PDF),
Power cord,
Dimensions: 41 x 37 x 15.5 cm,
Set weight: 12 kg,
Valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,

