



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DATEX-OHMEDA CARDIOCAP 5 patient monitor/cardiomonitor (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in the U.S.A.,  
Power supply: 230 V – no backup power supply,  
Frequency: 50/60 Hz,  
Power: 80 VA,  
Specifications:  
10.4" color LCD monitor,  
Built-in memory card slot,  
Clear and easy-to-use menu,  
RS 232 inputs, internet, REMOTE CONTROL, analog/digital output,  
Acoustic alarms,  
Spirometry analysis function  
Equipped with the following modules with cabling and sensors:  
ECG - 3-lead cable,  
SPO 2 - Clip sensor,  
NIBP - Cuff + airway,  
TEMP1 - (sensor not included - can be retrofitted with (additional charge)  
M-CAIO anesthetic gas analysis module with D-FEND water trap - CO2 measurement only,  
Menu language selection: ENG, GER, FR, ND,  
Includes:  
Power cord,  
NIBP cuff with airway,  
3-lead ECG sensor,  
SPO2 clip sensor,  
Dimensions: 33 x 22 x 30 cm,  
Weight: 11.2 kg,  
Updated inspection and ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional charge to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

