



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Microwave diathermy RADARTHERM S-60-D/HNO (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Technical specifications:
Timer switch: 0-30 minutes,
Smooth DOSIS adjustment: 0-60,
Continuous emission,
F=2450 +/-50 MHz,
2 applicators on movable arms ($\varnothing = 8$ cm) - applicator diameter,
Includes:
Power cord,
Overall dimensions: 25 x 25 x 110 cm,
Weight: 35 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

