



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Shortwave diathermy RADARTHERM S-60-D/HNO (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Technical specifications:  
Timer switch: 0-30 minutes,  
Smooth DOSIS adjustment: 0-60,  
Continuous emission,  
F=2450 +/-50 MHz,  
2 applicators on movable arms (Ø = 8 cm) - applicator diameter,  
Includes:  
Power cord,  
Overall dimensions: 25 x 25 x 110 cm,  
Weight: 35 kg,  
Updated inspection and ready to use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty: 6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

