



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RECOMED R12-200P microwave diathermy (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 700 VA,
Features:
Pulse and continuous operation,
Easy to use,
Processor-controlled,
Program memory,
High safety standard,
For all high-frequency and heat therapy applications,
Specifications:
F=2450 +/- 50 MHz,
Treatment duration: From 1 to 30 minutes,
Pulse operation: From 1 to 50 pulses,
Power regulation from 5 to 200 W,
Includes:
Power cord,
2 applicators,
Dimensions:
Overall: 40 x 30 x 79 cm,
Applicator no. 1: 44 x 12.5 cm,
Applicator No. 2: $\varnothing$ (diameter 6 cm),
Weight: 44.5 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

