



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDELA BASIC Surgical Suction Pump (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Switzerland,

Power supply: 230 V,

Frequency: 50/60 Hz,

The MEDELA BASIC Surgical Suction Unit

is a small surgical suction unit for suctioning secretions and fluids during procedures. It is primarily used in treatment rooms and hospital emergency departments (EDs).

Specifications:

Generated negative pressure: -90 kPa/-675 mmHg,

Max. capacity 50 l/min.,

Smooth suction power adjustment,

Built-in pressure gauge,

Equipped with a mounting rail for jars and additional accessories,

Includes:

Power cord,

Brand new autoclavable jar with a new disposable insert,

Condensate container,

Overall dimensions: 25 x 30 x 25 cm,

Weight: 25 kg,

Updated inspection and ready to use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

