



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Surgical mammal CHEIRON ADELA (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the Czech Republic,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 110 VA,
Surgical suction unit
CHEIRON ADELA

is a small surgical suction unit for suctioning secretions and fluids during procedures. It is primarily used in treatment rooms and hospital emergency departments (EDs).

Specifications:

Generated negative pressure: -100 kPa,
Smooth suction power adjustment,
Built-in pressure gauge,

Includes:

Power cord,
Brand new autoclavable jar with a new disposable insert,
Condensate container,
Overall dimensions: 38 x 27 x 24 cm,
Weight: 20 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

