



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE PANDA IRES Neonatal Incubator (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in the U.S.A. (2017),

Power supply: 230 V,

Frequency: 50/60 Hz,

GE Panda iRes

is a warming incubator for infants. Thanks to alarm muting, an easy-to-read control panel, and many other features, the GE Panda iRes ensures constant monitoring of the newborn's temperature.

Specifications:

Temperature control range: 34°C - 37°C,

Light intensity adjustment using a knob,

Heating power adjustment range: 0 - 100% (adjustable in 5% increments),

Built-in spotlight,

6.2" color display,

Tilt on both sides within a range of 12°,

Double-sided pull-out accessory drawers,

Electric adjustment within a range of 18 cm,

Built-in patient sensor connector – sensor not included (can be retrofitted for an additional fee),

Menu language selection: NOR, LIT, KOR, JP, IT, HUN, GREK, GER, FR, FI, PL, DUN, CS, CHIN, TUR SWE, ESP, RUS, RUM, PT, ENG, DUN

Includes:

Power cord,

Mattress,

Dimensions:

Landing surface: 68 x 50 cm,

Total dimensions: 54 x 100 x 190/208 cm,

Maximum load capacity: 40 kg,

Weight: 100 kg,

Updated inspection and ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan

