



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS INTELLIVUE MP30 patient monitor/cardiometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
German production (2014),
Power supply: 230 V,
Battery backup: via two built-in 11.1 V / 7.8 Ah batteries – batteries not included (new batteries can be retrofitted for an additional fee),
Frequency: 50/60 Hz,
The IntelliVue MP30
is a patient monitor that monitors multiple vital signs.
The panel-mounted screen is made of 10-inch TFT LCD technology and operates at SVGA resolution.
The touchscreen and navigation joystick allow for direct access to many functions using simple commands.
The comfortable handle and durable housing make it easy to carry.
Measurement trends can be presented graphically, combining up to three curves on a single graph, allowing for easy identification of changes.
patient's physiological condition,
Rapidly changing parameter trends can be updated per beat and up to four high-resolution trend segments can be displayed,
IntelliVue technology - access to information from various systems in the hospital network - pharmacy, laboratory, ECG management,
radiology PACS, server-based hospital systems, and others - and deliver it to the patient monitor,
No separate hard drive or fan - IntelliVue monitors are small, quiet, stable, and start quickly,
The navigation knob makes it easy to enter information and navigate the screen,
Acoustic and visual alarms,
LAN,
Patient monitor equipped with a PHILIPS IntelliVue X2 measurement module:
ECG /RESP - 3-lead cable,
SPO2
- Silicone sensor,
NBP - Cuff with air tube,
TEMP - Temperature connector socket,
(sensor not included, can be (can be retrofitted for an additional fee),
Equipment mounted on a mobile stand with a built-in accessory basket,
Menu in: Polish,
Includes:
User's manual in English (PDF),
Power cord,
3-lead ECG cable,
Silicone SPO2 sensor,
NIBP cuff with airway,
Overall dimensions: 60 x 140 cm,
Weight: 14 kg,
Cat 04,
Valid inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

