



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ASTAR VERSA equipment table (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Poland (2013),
Power supply: 230 V,
Frequency: 50/60 Hz,
ASTAR

is a Versa equipment table.

It stands out from the competition thanks to its intelligently designed shelves and spacious worktops that facilitate the organization of medical equipment. Whether you're using it in a small office or a large medical center, its versatility allows for efficient workspace management.

The table sits on a mobile base with 4 braking casters.

Dimensions:

Top shelf: 38 x 26 cm,

Bottom shelf: 36 x 25 cm,

Total dimensions: 65 x 50 x 93 cm,

Weight: 30 kg.

It has a current inspection and is ready to use.

A Technical Passport (Service Report) is issued, valid for 12 months.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

