



Surgical diathermy ERBE ICC 50 INT (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 50 W,
- **Technical specifications:**
 - Digital displays for cutting and coagulation power,
 - Foiled, easy-to-use control panel,
- **It has the following functions:**
 - Cut (CUT) - 50 W at 500 Ω:
 - EFFECT 1 (low coagulation effect),
 - EFFECT 2 (high coagulation effect),
 - Coagulation (COAG) - 50 W at 200 Ω:
 - Forced (FORCED),
 - Soft (SOFT),
 - Bipolar (BIPOLAR),
- **Included in the set:**
 - Monopolar electrode cable with ERGO PEN BOWA holder - 2 buttons, 2.4 mm shaft (**electrodes not included, can be retrofitted for an additional fee**),
 - BOWA clip cable for disposable neutral electrodes (**electrodes not included, can be retrofitted for an additional fee**),
 - Bipolar cable for BOWA tweezers (**can be retrofitted with bipolar tweezers and tweezers for an additional fee**),
 - User manual in English (PDF),
 - Power cord,
 - Wired foot switch,
- Dimensions: 27 x 26.5 x 10.5 cm,
- Weight: 4.5 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

