



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Surgical diathermy ERBE ICC-200 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 200 W,  
Smooth power regulation,  
Digital power setting displays,  
Laminate control panel,  
Control panel keypad in English,  
Features the following functions:  
ENDO-CUT,  
Auto-cutting - max. 200 W,  
Auto-coagulation - max. 120 W,  
Available coagulation modes: SOFT, FORCED, BIPOLAR,  
Included:  
ERBE foot switch,  
ERBE neutral electrode clip cable - electrodes not included, can be retrofitted for an additional fee,  
ERBE monopolar electrode holder (two buttons, 4 mm shaft),  
Bipolar cable (can be retrofitted with BOWA forceps and bipolar tweezers for an additional fee),  
Instruction manual in English (PDF),  
Power cord,  
Dimensions: 37 x 28 x 15 cm,  
Weight: 9.5 kg,  
Updated inspection and ready to use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee of up to 12 months for domestic and international markets. Financing options (Poland only): Installments, Leasing, Loan.

