



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SEIKO PRECISION VP-1200 videoprinter (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Japan,

Power supply: 230 V,

Frequency: 50/60 Hz,

The SEIKO PRECISION VP-1200 Video Printer

is a compact, lightweight, economical, and easy-to-connect and use thermal video printer.

Product features:

Connects to ultrasound scanners, electron microscopes, test, measurement, and analytical devices, as well as industrial and security monitoring systems.

256 shades of gray (8-bit).

Resolution: 848 x 490.

Print time: 2.9 – 3.5 s.

Print size:

Normal: 98 x 74 mm.

Wide: 104 x 77 mm.

The device automatically selects between image standards. NTSC and PAL,

Economical printing option provides up to 265 sheets per roll,

Easy access to contrast and brightness knobs and start, print, copy, and paper feed buttons on the front control panel,

Features BNC sockets for VIDEO IN-OUT and a stereo remote control jack,

Ability to change settings, including: Print size, direction, and length, sharpness, memory, impedance,

Choice of three paper types and three gamma types for high-contrast paper,

Includes:

Power cord,

Dimensions: 15.5 x 31 x 10 cm,

Weight: 3.4 kg,

Updated inspection and ready to use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

