



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SMITHS MEDICAL BCI 3301 pulse oximeter (Reconditioned)

Product temporarily unavailable
Used, in good condition,
Made in the USA,
Battery powered: 3 x R14,
Durable and lightweight construction,
Ergonomically designed to fit comfortably in the hand,
Bright and easy-to-read display,
Quantitative pulse strength bar on the screen,
Fast, reliable SpO2 and pulse rate measurements,
Stores data for up to 99 different patients,
Simple on/off control with automatic shutdown to save battery life,
Intended for patients of all ages, from newborns to adults,
Includes:
Clip sensor with extension cable,
3 x R14 batteries,
Dimensions: 8.5 x 17.5 x 4 cm,
Weight: 0.4 kg,
Valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
3 months for the market international,
Financing options (Poland only): Installments, Leasing, Loan,

