



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CV-165 video processor with CLE-165 light source (set) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Japan,
OLYMPUS CV-165 video processor:
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 150 VA,
Automatic white balance control,
Features "FREEZE" function,
Adjustable color saturation: red and blue,
Sealed control panel,
Built-in memory card reader,
Features the following inputs and outputs:
PRINTER REMOTE,
PC REMOTE,
MONITOR REMOTE,
RGB OUT,
Y/C OUT,
COMP. OUTPUT,
MONITOR OUTPUT,
DIGITAL OUTPUT,
LIGHT SOURCE,
LIGHT CONTROL
OLYMPUS CLE-165 Light Source:
Power Supply: 230 V,
Frequency: 50/60 Hz,
Power: 350 VA,
Suitable for flexible endoscopes,
Main and emergency lighting: 2 x 15 V/150 W halogen bulbs,
Easy bulb access from the side of the device,
Laminated control panel,
Automatic or manual image brightness adjustment,
Built-in air pump with electronic power control,
Included:
Power cord (2 pcs.),
Dimensions
:
OLYMPUS CV-165 Video Processor: 50 x 42 x 9.5 cm,
OLYMPUS CLE-165 Light Source: 49 x 37 x 14 cm cm,
Weight: 23.5 kg,
It has a current inspection and is ready to operate.
A Technical Passport (Service Report) is issued, valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
The warranty can be extended for an additional fee to 12 months for the domestic and international markets.
The warranty period does not cover light bulbs and other light sources installed in the device.
Financing options (Poland only): Installments, Leasing, Loan.

