



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE ERBOKRYO AE Cryosurgery Apparatus (Reconditioned)

Product temporarily unavailable

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 40 W,

Application:

Ophthalmology,

Dermatology,

Otolaryngology,

Proctology,

Veterinary Surgery,

Plastic surgery,

Gynecology,

Thoracic surgery,

Specifications:

Temperature setting range: -20 °C to -70 °C,

Defrost temperature: From 8 °C to 15 °C,

Cryogenic gas:

Gas type: Nitrous oxide (N₂O) or carbon dioxide (CO₂),

Operating pressure: 40 - 60 bar,

Gas consumption "freezing": 36 g/min,

Amount of gas produced: 40 - 60 l/min,

Instrumentation placed on ERBE mobile table:

Mobile 4-wheel (2 braking wheels),

Pull-out drawer for accessories: 48 x 32 x 8 cm,

Built-in cylinder mount,

Includes:

ERBE 20716-087 Ø 3.0 mm probe,

Dedicated ERBE mobile table,

Hose for supply cylinder,

Power cord,

Dimensions:

Appliances: 47 x 35 x 10 cm,

Overall: 64 x 50 x 90 cm,

Weight:

Devices: 10 KG,

Total: 25 KG,

It has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

